

## UHPLC separation of tocopherols

Method: VFD0092J  
HPLC: RP UHPLC

Column: Eurospher II 100-2 C18 A, 100 x 2 mm ID Order No. 10BE184E2F

Phase: Eurospher II 100-2 C18 A

Conditions:

Euent:	A: Water
	B: ACN
Gradient:	0-4 min 80%-95% B
	4-6 min 95% B
	6.02 min 80% B, 3 min hold
Flow rate:	0.6 ml/min
Temperature:	30 °C
Volume:	2 µL

Detection: 215 nm

Substances: Alpha-, gamma-, delta-tocopherol

Keywords: tocopherols

## Chromatogram:

1 Imp. A  
2 delta Tocopherol  
3 gamma Tocopherol  
4 alpha Tocopherol

