




Read the corresponding technical documentation for handling and safety reasons.

## Running-in information



After installing a new deuterium lamp into your KNAUER detector, ensure to allow the lamp a running-in time of approximately 24 hours.

For further details on installing the lamp, consult the technical documentation of the detector.

Lamp replacement, V6757, Version 2.0



**KNAUER**

**KNAUER**

**Wissenschaftliche Geräte GmbH**

**Hegauer Weg 38**

**14163 Berlin, Germany**

**Phone: +49 30 809727-0**

**Fax: +49 30 8015010**

**E-Mail: [info@knauer.net](mailto:info@knauer.net)**

**Internet: [www.knauer.net](http://www.knauer.net)**